

PRODUCT SPECIFICATION

SAFETYWARE LatexTouch™ 9" Powder Free Latex Examination Gloves

- Medical grade of AQL 2.5, which is suitable for healthcare and dental applications.
- Made of elastic and soft latex for superior dexterity.
- Polymer Coated for smooth donning and doffing experience.
- Textured grip for secure handling and precision.
- Enhanced tactile sensitivity for delicate tasks.

APPLICATIONS

- Laboratory
- Cleaning Work
- General Industry

STORAGE

- Store in a cool, dry place. Avoid excessive heat.
- Open box should be shielded from exposure to direct sun or fluorescent lighting.

TECHNICAL DATA

MATERIAL	Powder free latex	
SURFACE TEXTURE	Beaded Cuff, Ambidextrous and Palm Textured.	
COLOR	Natural white	
SIZE	S, M, L, XL	
THICKNESS (SINGLE WALL)	Finger : 0.15 ± 0.02 mm Palm : 0.13 ± 0.02 mm Cuff : 0.10 ± 0.02 mm	
	Before Aging	After Aging
TENSILE STRENGTH (Mpa)	min 18	min 14
ELONGATION AT BREAK (%)	650	500
LENGTH	Min 240 mm	
STANDARD	ASTM D 6124, ASTM D 5712, ASTM D 3578, EN 455 CE 2777	
PACKAGING	100pcs /box, 10 boxes /carton	
SHELF LIFE	3 Years	

Item Code: LEX55



ITEM CODE

SIZE	ITEM CODE	WEIGHT (± 5.0 g)	PALM WIDTH (± 5.0 mm)
S	LEX55-SM	5.5	85
M	LEX55-MD	6.0	95
L	LEX55-LG	6.5	105
XL	LEX55-XL	7.0	115

SAFETYWARE SOLUTIONS SDN. BHD. (99556-U)
Plot 237, Lengkok Perindustrian Bukit Minyak 3,
Bukit Minyak Industrial Estate,
14100 Simpang Ampat, Penang, Malaysia.

Tel : +604-5023 882 | Fax : +604-5081 882
www.safetyware.com | info@safetyware.com

It is the product users' responsibility to conduct all the appropriate evaluations & testing to determine the suitability of SAFETYWARE products for use in a particular environment or purpose. SAFETYWARE disclaims all warranties other than as expressly provided.

Date created: 14/11/2023

Revision: 01

Date revised: 04/06/2026