

# PRODUCT SPECIFICATION

## SAFETYWARE NiteTouch NEX35 9" Powder Free Light Duty Nitrile Examination Gloves

- Crafted from high-quality nitrile material.
- Designed for light-duty tasks, these gloves are ideal for non-sterile examination procedures, hygiene-sensitive environments, and short-duration work that still requires dependable protection.
- Offer enhanced resistance to punctures, tears, and chemicals compared to traditional latex or vinyl gloves without latex allergens.
- Featuring textured fingertips for precision control and a powder-free formulation to prevent skin irritation and contamination,

### APPLICATIONS

- Healthcare & Dental
- Laboratory & Research
- Cleaning Work
- General Industry
- Food Processing & Service

### TECHNICAL DATA

|                        |                                |             |
|------------------------|--------------------------------|-------------|
| MATERIAL               | Nitrile (Latex-Free)           |             |
| SURFACE TEXTURE        | Beaded Cuff, Finger Textured   |             |
| COLOR                  | Blue                           |             |
| SIZE                   | S, M, L, XL                    |             |
| WEIGHT                 | 3.5 ± 0.3 g                    |             |
| LENGTH                 | 240 ± 5 mm                     |             |
|                        | Before Aging                   | After Aging |
| TENSILE STRENGTH (Mpa) | Min 14                         | Min 14      |
| ELONGATION(%) (min)    | 500                            | 400         |
| STANDARD               | ASTM D6319, EN 455             |             |
| PACKAGING              | 100 pcs /box, 10 boxes /carton |             |
| SHELF LIFE             | 5 Years                        |             |

Item Code: NEX35



### ITEM CODE

| SIZE | ITEM CODE | PALM WIDTH (± 5.0 mm) |
|------|-----------|-----------------------|
| S    | NEX35-SM  | 85                    |
| M    | NEX35-MD  | 95                    |
| L    | NEX35-LG  | 105                   |
| XL   | NEX35-XL  | 115                   |

### STORAGE

- Store in a cool, dry place. Avoid excessive heat.
- Open box should be shielded from exposure to direct sun or fluorescent lighting.